

ULTRASOUND WED-M1 vet

Digital Veterinary Scanner

Digital WED-M1 vet, is the new contribution by WELLD for veterinary in 2008. It has the transducer, display and power supply integrated in one compact sturdy housing. The total unit weights 500g (including the battery)only. After use, it can be simply be cleaned under running tap water.



CE

 **WELLD**
ULTRASOUND SCANNER

WED-M1 vet

Digital veterinary Scanner

System Feature

- Digital technology
- 3.5 inch TFT color LCD
- Washable
- Waterproof(Main Unit--XP4,Probe--XP7)
- Simple to operate
- All in one, no cable
- 3hrs continuous operation
- 2 years warranty

Physical Characteristic

- Size: 171×93×62 mm
- Weight: 500g including the battery

Standard Accessories

- AC-adapter (110/220V)
- two rechargeable battery
- One battery charging cable
- One carrying bag
- One bottle of ultra gel
- One hand strap



Specifications

Monitor : digital 3.5 inch TFT color LCD

Probe : Multi-frequency mechanical sector probe

Scan Frequency : 2.5MHz, 3.5MHz, 5.0MHz

Display Depth : MAX--160mm

Scan Angle : 78 degree

Gray Shade : 256 levels

Measurement : Distance, Area, Volume, Heart Rate

OB Measurment : Swine-HL, Equine, Bovine-BL,SL,HL, Sheep-USD

Cat-HD,BD, Dog-BD, GSD, CRL,HD, BD

Permanent Storage : 128 frames

AC-adapter : 110V~240V, 50/60Hz

Battery : (1550mAh) Internal rechargeable battery

More than 3hrs for continuous operation

less than 3hrs for recharging



Shenzhen WELL.D Medical Electronics Co.,Ltd.

<http://www.welld.net> <http://www.welld.com.cn> Email: export@welld.com.cn

ADD:12/F BLock North, New Energy Bldg., Nanhai Ave., Nanshan District, ShenZhen 518054,Guangdong, China
Tel:+86-755-26072063/805/821 Fax:+86-755-26073350